



**Animal Clinic
Northview, Inc.**

36400 Center Ridge Rd.
N. Ridgeville, OH 44039
repro@northviewvet.com
Phone: (440) 327-8282
Fax: (440) 327-8845

Sample Preparation

Draw 3 mL/cc of blood in a plain red top tube (no serum separator). Centrifuge, draw off serum, and send in another plain Red top tube also a Green top Li-Heparin LH/1.3 tube with an ice pack overnight. Use overnight courier of choice , i.e. Federal Express,UPS, US Postal Service.

Please send an unstained vaginal smear with all first progesterone samples.

Progesterone Sample Submission Form

Name of Dog: _____

Breed of Dog: _____

DOB: _____ Color of Dog: _____

Owner: _____

Address: _____

City: _____ State: _____ Zip Code: _____

Phone: _____

There will be a one time interpretation fee charged if the breeding is not performed here.

Credit Card Number: _____

Expiration Date: _____ V-Code: _____

Check Enclosed: _____

Signature of Card Holder: _____

Date of Sample: _____

Type of Semen: Fresh Fresh Chilled Frozen

Insemination Type: Vaginal Surgical

Transcervical Natural

Person To Be Notified: _____

Phone: _____